

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Dog GFAP(Glial Fibrillary Acidic Protein) ELISA Kit EBT-ELK9719

Artikelname	Dog GFAP(Glial Fibrillary Acidic Protein) ELISA Kit
Artikelnummer	EBT-ELK9719
Hersteller Artikelnummer	ELK9719
Alternativnummer	EBT-ELK9719-96, EBT-ELK9719-48, EBT-ELK9719-96X5
Hersteller	ELK Biotechnology
Kategorie	Kits/Assays
Spezies Reaktivität	Canine
Konzentration	10 ng/mL
Detektionsbereich	0.16-10 ng/mL
Sensitivitaet	0.055 ng/mL
Proben	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Anwendungsbeschreibung	Assay Type: Sandwich. Assay length: 3.5h. Research Area: Infection immunity, Neuro science,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Dog GFAP. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Dog GFAP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Dog GFAP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Dog GFAP in the samples is then determined by comparing the OD of the samples to the standard curve
------------------------	---