

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human SUMO-conjugating enzyme UBC9 (UBE2I/UBC9/UBCE9) ELISA Kit ASC-KTE60101

Artikelname	Human SUMO-conjugating enzyme UBC9 (UBE2I/UBC9/UBCE9) ELISA Kit
Artikelnummer	ASC-KTE60101
Hersteller Artikelnummer	KTE60101
Alternativnummer	ASC-KTE60101-48, ASC-KTE60101-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human SUMO-conjugating enzyme UBC9 (UBE2I/UBC9/UBCE9) ELISA Kit employs a two-site sandwich ELISA to quantitate UBE2I....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	UBE2I
NCBI	7329
UniProt	P63279
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Human SUMO-conjugating enzyme UBC9 (UBE2I/UBC9/UBCE9) ELISA Kit employs a two-site sandwich ELISA to quantitate UBE2I in samples. An antibody specific for UBE2I has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any UBE2I present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for UBE2I is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of UBE2I bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	---