

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Human TNF receptor-associated factor 3 (TRAF3) ELISA Kit ASC-KTE60140

Artikelname	Human TNF receptor-associated factor 3 (TRAF3) ELISA Kit
Artikelnummer	ASC-KTE60140
Hersteller Artikelnummer	KTE60140
Alternativnummer	ASC-KTE60140-48, ASC-KTE60140-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human TNF receptor-associated factor 3 (TRAF3) ELISA Kit employs a two-site sandwich ELISA to quantitate TRAF3....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TRAF3
NCBI	<a href="#">7187</a>
UniProt	<a href="#">Q13114</a>
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

This Human TNF receptor-associated factor 3 (TRAF3) ELISA Kit employs a two-site sandwich ELISA to quantitate TRAF3 in samples. An antibody specific for TRAF3 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TRAF3 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TRAF3 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TRAF3 bound in the initial step. The color development is stopped and the intensity of the color is measured.