

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Human Taste receptor type 2 member 40 (TAS2R40) ELISA Kit ASC-KTE60499

Artikelname	Human Taste receptor type 2 member 40 (TAS2R40) ELISA Kit
Artikelnummer	ASC-KTE60499
Hersteller Artikelnummer	KTE60499
Alternativnummer	ASC-KTE60499-48, ASC-KTE60499-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Taste receptor type 2 member 40 (TAS2R40) ELISA Kit employs a two-site sandwich ELISA to quantitate TAS2R40....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TAS2R40
NCBI	<a href="#">259286</a>
UniProt	<a href="#">P59535</a>
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

This Human Taste receptor type 2 member 40 (TAS2R40) ELISA Kit employs a two-site sandwich ELISA to quantitate TAS2R40 in samples. An antibody specific for TAS2R40 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TAS2R40 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TAS2R40 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TAS2R40 bound in the initial step. The color development is stopped and the intensity of the color is measured.