

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Human Mothers Against Decapentaplegic Homolog 3 (SMAD3) ELISA Kit ASC-KTE60574

Artikelname	Human Mothers Against Decapentaplegic Homolog 3 (SMAD3) ELISA Kit
Artikelnummer	ASC-KTE60574
Hersteller Artikelnummer	KTE60574
Alternativnummer	ASC-KTE60574-48, ASC-KTE60574-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Mothers Against Decapentaplegic Homolog 3 (SMAD3) ELISA Kit employs a two-site sandwich ELISA to quantitate SMAD3....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	SMAD3
NCBI	<a href="#">4088</a>
UniProt	<a href="#">P84022</a>
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

This Human Mothers Against Decapentaplegic Homolog 3 (SMAD3) ELISA Kit employs a two-site sandwich ELISA to quantitate SMAD3 in samples. An antibody specific for SMAD3 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any SMAD3 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SMAD3 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SMAD3 bound in the initial step. The color development is stopped and the intensity of the color is measured.