

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Myosin-IIIb (MYO3B) ELISA Kit ASC-KTE61417

| | |
|--------------------------|--|
| Artikelname | Human Myosin-IIIb (MYO3B) ELISA Kit |
| Artikelnummer | ASC-KTE61417 |
| Hersteller Artikelnummer | KTE61417 |
| Alternativnummer | ASC-KTE61417-48, ASC-KTE61417-96 |
| Hersteller | Abbkine Scientific |
| Kategorie | Kits/Assays |
| Applikation | ELISA |
| Spezies Reaktivität | Human |
| Produktbeschreibung | This Human Myosin-IIIb (MYO3B) ELISA Kit employs a two-site sandwich ELISA to quantitate MYO3B.... |
| Detektionsbereich | Please inquire |
| Sensitivitaet | Please inquire |
| Tag | MYO3B |
| NCBI | 140469 |
| UniProt | Q8WXR4 |
| Proben | Cell culture supernatantsSerumPlasmaOther biological fluids |

Anwendungsbeschreibung

This Human Myosin-IIIb (MYO3B) ELISA Kit employs a two-site sandwich ELISA to quantitate MYO3B in samples. An antibody specific for MYO3B has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any MYO3B present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for MYO3B is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of MYO3B bound in the initial step. The color development is stopped and the intensity of the color is measured.