

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Human MAP kinase-activated protein kinase 5 (MAPKAPK5) ELISA Kit ASC-KTE61728

Artikelname	Human MAP kinase-activated protein kinase 5 (MAPKAPK5) ELISA Kit
Artikelnummer	ASC-KTE61728
Hersteller Artikelnummer	KTE61728
Alternativnummer	ASC-KTE61728-48, ASC-KTE61728-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human MAP kinase-activated protein kinase 5 (MAPKAPK5) ELISA Kit employs a two-site sandwich ELISA to quantitate MAPKAPK5....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	MAPKAPK5
NCBI	<a href="#">8550</a>
UniProt	<a href="#">Q8IW41</a>
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

This Human MAP kinase-activated protein kinase 5 (MAPKAPK5) ELISA Kit employs a two-site sandwich ELISA to quantitate MAPKAPK5 in samples. An antibody specific for MAPKAPK5 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any MAPKAPK5 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for MAPKAPK5 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of MAPKAPK5 bound in the initial step. The color development is stopped and the intensity of the color is measured.