

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse Tenascin-R (TNR) ELISA Kit ASC-KTE70090

Artikelname	Mouse Tenascin-R (TNR) ELISA Kit
Artikelnummer	ASC-KTE70090
Hersteller Artikelnummer	KTE70090
Alternativnummer	ASC-KTE70090-48, ASC-KTE70090-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse Tenascin-R (TNR) ELISA Kit employs a two-site sandwich ELISA to quantitate TNR....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TNR
NCBI	21960
UniProt	Q8BYI9
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Mouse Tenascin-R (TNR) ELISA Kit employs a two-site sandwich ELISA to quantitate TNR in samples. An antibody specific for TNR has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNR present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNR is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNR bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--