

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse Meckelin (TMEM67) ELISA Kit ASC-KTE70164

| | |
|--------------------------|---|
| Artikelname | Mouse Meckelin (TMEM67) ELISA Kit |
| Artikelnummer | ASC-KTE70164 |
| Hersteller Artikelnummer | KTE70164 |
| Alternativnummer | ASC-KTE70164-48, ASC-KTE70164-96 |
| Hersteller | Abbkine Scientific |
| Kategorie | Kits/Assays |
| Applikation | ELISA |
| Spezies Reaktivität | Mouse |
| Produktbeschreibung | This Mouse Meckelin (TMEM67) ELISA Kit employs a two-site sandwich ELISA to quantitate TMEM67.... |
| Detektionsbereich | Please inquire |
| Sensitivitaet | Please inquire |
| Tag | TMEM67 |
| NCBI | 329795 |
| UniProt | Q8BR76 |
| Proben | Cell culture supernatantsSerumPlasmaOther biological fluids |

Anwendungsbeschreibung

This Mouse Meckelin (TMEM67) ELISA Kit employs a two-site sandwich ELISA to quantitate TMEM67 in samples. An antibody specific for TMEM67 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TMEM67 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TMEM67 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TMEM67 bound in the initial step. The color development is stopped and the intensity of the color is measured.