

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse Tight junction protein ZO-2 (TJP2) ELISA Kit ASC-KTE70219

Artikelname	Mouse Tight junction protein ZO-2 (TJP2) ELISA Kit
Artikelnummer	ASC-KTE70219
Hersteller Artikelnummer	KTE70219
Alternativnummer	ASC-KTE70219-48, ASC-KTE70219-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse Tight junction protein ZO-2 (TJP2) ELISA Kit employs a two-site sandwich ELISA to quantitate TJP2....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TJP2
NCBI	21873
UniProt	Q9Z0U1
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

This Mouse Tight junction protein ZO-2 (TJP2) ELISA Kit employs a two-site sandwich ELISA to quantitate TJP2 in samples. An antibody specific for TJP2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TJP2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TJP2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TJP2 bound in the initial step. The color development is stopped and the intensity of the color is measured.