

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Mouse TIMELESS-interacting protein (TIPIN) ELISA Kit ASC-KTE70222

Artikelname	Mouse TIMELESS-interacting protein (TIPIN) ELISA Kit
Artikelnummer	ASC-KTE70222
Hersteller Artikelnummer	KTE70222
Alternativnummer	ASC-KTE70222-48, ASC-KTE70222-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse TIMELESS-interacting protein (TIPIN) ELISA Kit employs a two-site sandwich ELISA to quantitate TIPIN....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TIPIN
NCBI	<a href="#">66131</a>
UniProt	<a href="#">Q91WA1</a>
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

This Mouse TIMELESS-interacting protein (TIPIN) ELISA Kit employs a two-site sandwich ELISA to quantitate TIPIN in samples. An antibody specific for TIPIN has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TIPIN present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TIPIN is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TIPIN bound in the initial step. The color development is stopped and the intensity of the color is measured.