

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse Testican-2 (SPOCK2) ELISA Kit ASC-KTE70236

Artikelname	Mouse Testican-2 (SPOCK2) ELISA Kit
Artikelnummer	ASC-KTE70236
Hersteller Artikelnummer	KTE70236
Alternativnummer	ASC-KTE70236-48, ASC-KTE70236-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse Testican-2 (SPOCK2) ELISA Kit employs a two-site sandwich ELISA to quantitate SPOCK2....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	SPOCK2
NCBI	94214
UniProt	Q9ER58
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Mouse Testican-2 (SPOCK2) ELISA Kit employs a two-site sandwich ELISA to quantitate SPOCK2 in samples. An antibody specific for SPOCK2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any SPOCK2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SPOCK2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SPOCK2 bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--