

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ExKine(TM) Membrane and Cytoplasmic Protein Extraction kit ASC-KTP3005

Artikelname	ExKine(TM) Membrane and Cytoplasmic Protein Extraction kit
Artikelnummer	ASC-KTP3005
Hersteller Artikelnummer	KTP3005
Alternativnummer	ASC-KTP3005-50, ASC-KTP3005-200
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Produktbeschreibung	Abbkine ExKine(TM) Membrane and Cytoplasmic Protein Extraction Kit enable stepwise separation and preparation of cytoplasmic and membrane protein....
Tag	Cytoplasmic
Anwendungsbeschreibung	Abbkine ExKine(TM) Membrane and Cytoplasmic Protein Extraction Kit enable stepwise separation and preparation of non-denatured, active cytoplasmic and membrane protein from mammalian cultured cells or tissues, which can be used directly in a variety of proteomics application. The cells are first broken by homogenization to allow the release of cytosolic and membrane proteins, after which membrane proteins are obtained by high speed centrifugation. Extraction efficiencies and yields will vary depending on cell type as well as the number of times the integral membrane protein spans the lipid bilayer.