

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit RDD-RDR-NT-PROANP-MU

Artikelname	Mouse N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit
Artikelnummer	RDD-RDR-NT-PROANP-MU
Hersteller Artikelnummer	RDR-NT-ProANP-Mu
Alternativnummer	RDD-RDR-NT-PROANP-MU-96,RDD-RDR-NT-PROANP-MU-48
Hersteller	reddot Biotech
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit is a precoated, Sandwich ELISA assay designed to detect and quantify the target analyte in serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other...
Detektionsbereich	78.1-5000pg/mL
Sensitivitaet	29pg/mL
Proben	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids
Target-Kategorie	Mouse N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP)
Anwendungsbeschreibung	Assay Length: 3 hours. Experimental Method: Sandwich