

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Equine Procollagen II N-Terminal Propeptide (PIINP) ELISA Kit RDD-RDR-PIINP-EQ

| | |
|--------------------------|---|
| Artikelname | Equine Procollagen II N-Terminal Propeptide (PIINP) ELISA Kit |
| Artikelnummer | RDD-RDR-PIINP-EQ |
| Hersteller Artikelnummer | RDR-PIINP-Eq |
| Alternativnummer | RDD-RDR-PIINP-EQ-96,RDD-RDR-PIINP-EQ-48 |
| Hersteller | reddot Biotech |
| Kategorie | Kits/Assays |
| Applikation | ELISA |
| Spezies Reaktivität | Equine |
| Produktbeschreibung | This Equine Procollagen II N-Terminal Propeptide (PIINP) ELISA Kit is a precoated, Sandwich ELISA assay designed to detect and quantify the target analyte in serum, plasma and other biological fluids . This product is for Research Use Only.... |
| Detektionsbereich | 3.12-200ng/mL |
| Sensitivitaet | 1.34ng/mL |
| Proben | serum, plasma and other biological fluids |
| Target-Kategorie | Equine Procollagen II N-Terminal Propeptide (PIINP) |
| Anwendungsbeschreibung | Assay Length: 3 hours. Experimental Method: Sandwich |