

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

MAP Kinase p44/42 (MAPK3/MAPK1, Mitogen-activated Protein Kinase 1, MAP Kinase 1, MAPK 1, Mapk, ERT1, Extracellular Signal-regulated Kinase 2, ERK-2, MAP Kinase Isoform p42, p42-MAPK, Mitogen-activated Protein Kinase 2, MAP Kinase USB-M2360-02G-ML650

Artikelname	MAP Kinase p44/42 (MAPK3/MAPK1, Mitogen-activated Protein Kinase 1, MAP Kinase 1, MAPK 1, Mapk, ERT1, Extracellular Signal-regulated Kinase 2, ERK-2, MAP Kinase Isoform p42, p42-MAPK, Mitogen-activated Protein Kinase 2, MAP Kinase
Artikelnummer	USB-M2360-02G-ML650
Hersteller Artikelnummer	M2360-02G-ML650
Alternativnummer	USB-M2360-02G-ML650-100
Hersteller	US Biological
Wirt	Rabbit
Kategorie	Antikörper
Applikation	FLISA, IHC, IP, WB
Immunogen	Synthetic peptide corresponding to aa339-353 from the C-terminal segment of rat ERK1.
Produktbeschreibung	MaxLight(TM) 650 is a new Far-IR stable dye conjugate comparable to Alexa Fluor(TM)647, DyLight(TM)649, Cy5(TM) and offers better labeling efficiency, brighter imaging and increased immunodetection. Absorbance (655nm), Emission (676nm), Extinction Co...
Reinheit	Purified by immunoaffinity chromatography.

Formulierung	Supplied as a liquid in PBS, pH 7.2. No preservatives added. Labeled with MaxLight(TM)650.
--------------	--