

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Mouse Nitric Oxide Synthase 2, Inducible (NOS2/iNOS) Microsample Fast ELISA Kit BYT-ORB1806587

Artikelname	Mouse Nitric Oxide Synthase 2, Inducible (NOS2/iNOS) Microsample Fast ELISA Kit
Artikelnummer	BYT-ORB1806587
Hersteller Artikelnummer	orb1806587
Alternativnummer	BYT-ORB1806587-48,BYT-ORB1806587-96
Hersteller	Biorbyt
Kategorie	Kits/Assays
Spezies Reaktivität	Mouse
Produktbeschreibung	Mouse Nitric Oxide Synthase 2, Inducible (NOS2/iNOS) Microsample Fast ELISA Kit...
Detektionsbereich	0.32-20ng/mL
Sensitivitaet	0.19 ng/mL
UniProt	<a href="#">P29477</a>
Formulierung	Ready to use
Proben	serum, plasma, Tissue homogenate and Other biological samples,Sample volume:25µL
Target-Kategorie	NOS2/iNOS

Anwendungsbeschreibung

Application Notes: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse NOS2/iNOS. Samples (or Standards) are added to the micro ELISA plate wells and combined with the specific antibody and biotinylated detection antibody specific for Mouse NOS2/iNOS. Then Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Mouse NOS2/iNOS, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm. The OD value is proportional to the concentration of Mouse NOS2/iNOS. You can calculate the concentration of Mouse NOS2/iNOS in the samples by comparing the OD of the samples to the standard curve