

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### FLICA 660 Poly Caspase Assay Kit BYT-ORB1409918

Artikelname	FLICA 660 Poly Caspase Assay Kit
Artikelnummer	BYT-ORB1409918
Hersteller Artikelnummer	orb1409918
Alternativnummer	BYT-ORB1409918-25
Hersteller	Biorbyt
Kategorie	Kits/Assays
Produktbeschreibung	This in vitro assay employs the far-red fluorescent inhibitor probe 660-VAD-FMK to label active caspase enzymes in living cells. Analyze the fluorescent signal using fluorescence microscopy or flow cytometry....
Proben	Cell culture
Anwendungsbeschreibung	Application Notes: 1. Prepare samples and controls., 2. Dilute 10X Cellular Wash Buffer 1:10 with diH2O., 3. Reconstitute FLICA with 50 µL DMSO., 4. Dilute FLICA 1:5 by adding 200 µL PBS., 5. Add diluted FLICA to each sample at 1:30-1:60 (e.g. spike at 1:30 by adding 10 µL to 290 µL cultured cells)., 6. Incubate approximately 1 hour., 7. Remove media and wash cells 3 times: add 1X Cellular Wash Buffer and spin cells., 8. If desired, label with additional stains, such as Hoechst, DAPI, or an antibody., 9. If desired, fix cells., 10. Analyze with a fluorescence microscope or flow cytometer. FLICA 660 excites at 660 nm and emits at 680-690 nm