

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Cell-mediated Cytotoxicity Assay: Basic Cytotoxicity Test BYT-ORB1409926

Artikelname	Cell-mediated Cytotoxicity Assay: Basic Cytotoxicity Test
Artikelnummer	BYT-ORB1409926
Hersteller Artikelnummer	orb1409926
Alternativnummer	BYT-ORB1409926-125
Hersteller	Biorbyt
Kategorie	Sonstiges
Produktbeschreibung	A dual-color assay for determining cell toxicity by flow cytometry. A green fluorescent cellular stain, CFSE, labels target cells and the red live/dead viability dye 7-AAD to identify the dead cells present in the cytotoxicity assay samples....
Proben	Cell culture
Anwendungsbeschreibung	Application Notes: 1. Prepare samples and controls, 2. Dilute Assay Buffer 1:10 with diH2O and filter sterilize., 3. Reconstitute CFSE with 200 µL DMSO and dilute 1:250 with 1X Assay Buffer., 4. Add 200 µL CFSE to target cells to stain them green and incubate 15 minutes at room temperature., 5. Wash target cells and adjust cell concentration., 6. Add effector cells to target cells., 7. Reconstitute 7-AAD with 260 µL DMSO and dilute 1:10 with 1X Assay Buffer., 8. Add 20 µL 7-AAD to label necrotic cells red and incubate 10 minutes on ice., 9. Run controls and analyze with a flow cytometer. Identify target cells by analyzing FL-1 (green) vs. FL-3 (red). CFSE excites at 492 nm and emits at 520-540 nm, 7-AAD excites at 546 nm and emits at 647 nm