

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Luciferase Antibody Peroxidase Conjugated, HRP, Goat, Polyclonal Preis auf Anfrage BYT-ORB344242

Artikelname	Luciferase Antibody Peroxidase Conjugated, HRP, Goat, Polyclonal Preis auf Anfrage
Artikelnummer	BYT-ORB344242
Hersteller Artikelnummer	orb344242
Alternativnummer	BYT-ORB344242-25
Hersteller	Biorbyt
Wirt	Goat
Kategorie	Antikörper
Applikation	ELISA, WB
Immunogen	Luciferase [Photinus pyralis (Firefly)]
Konjugation	HRP
Produktbeschreibung	Luciferase antibody (Peroxidase)...
Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Puffer	Preservative: 0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!. Stabilizer: 10 mg/mL Bovine Serum Albumin (rAlbumin) - Immunoglobulin and Protease free, Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Reinheit	<p>Anti-Luciferase Peroxidase Conjugated is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Goat Serum as well as purified and partially purified Luciferase [Photinus pyralis (Firefly)]. No reactivity is observed against Sea pansy (Renilla reniformis) luciferase.</p>
Formulierung	Liquid (sterile filtered)
Application Verdünnung	ELISA: 1:2,000 - 1:30,000, WB: 1:500 - 1:2,000
Anwendungsbeschreibung	<p>Application Notes: Anti-Luciferase Peroxidase conjugated antibody is produced from Photinus pyralis (Firefly) which produces a green light with a wavelength of 562 nm. Anti-Luciferase Peroxidase has been tested by ELISA and western blot. Expect a band at approximately 60.7 kDa in size corresponding to Luciferase by western blotting in the appropriate cell lysate or extract. Anti-Luciferase has been assayed against 1.0 ug of Luciferase [Photinus pyralis (Firefly)] in a standard capture ELISA using ABTS (2,2-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature. A working dilution of 1:50,000 to 1:200,000 of the reconstitution concentration is suggested for this product</p>