

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

BUB1 Antibody, IgG1, Clone: [14H5], Unconjugated, Mouse, Monoclonal Preis auf Anfrage BYT-ORB344427

Artikelname	BUB1 Antibody, IgG1, Clone: [14H5], Unconjugated, Mouse, Monoclonal Preis auf Anfrage
Artikelnummer	BYT-ORB344427
Hersteller Artikelnummer	orb344427
Alternativnummer	BYT-ORB344427-100
Hersteller	Biorbyt
Wirt	Mouse
Kategorie	Antikörper
Applikation	ELISA, IF, IP, WB
Spezies Reaktivität	Human
Immunogen	This protein A purified monoclonal antibody was produced by repeated immunizations with a recombinant protein corresponding to amino acid residues 281-419 of human BUB1 protein.
Konjugation	Unconjugated
Produktbeschreibung	BUB1 antibody...
Klonalität	Monoclonal
Konzentration	1.0 mg/mL
Klon-Bezeichnung	[14H5]

Isotyp	IgG1
NCBI	004327
UniProt	O43683
Puffer	Preservative: 0.01% (w/v) Sodium Azide. Stabilizer: None, Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	This Protein A purified antibody is directed against human BUB1 protein. The product was purified from tissue culture supernatant by chromatography. This antibody reacts with BUB1 from human cells. No reactivity is seen for homologues from Xenopus or Kangaroo Rat. Reactivity against homologues from other sources is not known.
Formulierung	Liquid (sterile filtered)
Application Verdünnung	ELISA: 1:5,000 - 1:20,000, IF: 1:500 - 1:2,000, IP: 1:200, WB: 1:500 - 1:2,000
Anwendungsbeschreibung	Application Notes: This protein A purified antibody has been tested for use in immunoprecipitation, immunofluorescence staining and western blot and is capable of detecting endogenous protein. Specific conditions for reactivity should be optimized by the end user. Expect a predominant band at ~ 160 kDa corresponding to full-length protein by western blotting in the appropriate cell lysate or extract. The use of HeLa whole cell lysates prepared using a RIPA buffer is recommended as a positive control. For IF microscopy use cells grown on cover slips fixed with 3.5% paraformaldehyde in PBS at pH 6.8 OR 100% methanol at -20 C. Permeabilize fixed cells with 0.5% Triton X-100