

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Pro Atrial Natriuretic Peptide (Pro-ANP) ELISA Kit BYT-ORB3146860

Artikelname	Human Pro Atrial Natriuretic Peptide (Pro-ANP) ELISA Kit
Artikelnummer	BYT-ORB3146860
Hersteller Artikelnummer	orb3146860
Alternativnummer	BYT-ORB3146860-48, BYT-ORB3146860-5, BYT-ORB3146860-50, BYT-ORB3146860-96
Hersteller	Biorbyt
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Pro Atrial Natriuretic Peptide (Pro-ANP) ELISA Kit employs a two-site sandwich ELISA to quantitate Pro-ANP....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

Application Notes: This Human Pro Atrial Natriuretic Peptide (Pro-ANP) ELISA Kit employs a two-site sandwich ELISA to quantitate Pro-ANP in samples. An antibody specific for Pro-ANP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any Pro-ANP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for Pro-ANP is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of Pro-ANP bound in the initial step. The color development is stopped and the intensity of the color is measured