

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Pyridoxal/pyridoxine vitamin B6-kinase (PDXK) ELISA Kit BYT-ORB3146981

Artikelname	Human Pyridoxal/pyridoxine vitamin B6-kinase (PDXK) ELISA Kit
Artikelnummer	BYT-ORB3146981
Hersteller Artikelnummer	orb3146981
Alternativnummer	BYT-ORB3146981-48, BYT-ORB3146981-5, BYT-ORB3146981-50, BYT-ORB3146981-96
Hersteller	Biorbyt
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Pyridoxal/pyridoxine vitamin B6-kinase (PDXK) ELISA Kit employs a two-site sandwich ELISA to quantitate PDXK....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

Application Notes: This Human Pyridoxal/pyridoxine vitamin B6-kinase (PDXK) ELISA Kit employs a two-site sandwich ELISA to quantitate PDXK in samples. An antibody specific for PDXK has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any PDXK present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PDXK is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of PDXK bound in the initial step. The color development is stopped and the intensity of the color is measured