

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit BYT-ORB3147762

Artikelname	Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit
Artikelnummer	BYT-ORB3147762
Hersteller Artikelnummer	orb3147762
Alternativnummer	BYT-ORB3147762-48, BYT-ORB3147762-5, BYT-ORB3147762-50, BYT-ORB3147762-96
Hersteller	Biorbyt
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit employs a two-site sandwich ELISA to quantitate TNIK....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung

Application Notes: This Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit employs a two-site sandwich ELISA to quantitate TNIK in samples. An antibody specific for TNIK has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTNIK present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNIK is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNIK bound in the initial step. The color development is stopped and the intensity of the color is measured