

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

General 2-Arachidonoylglycerol (2-AG) ELISA Kit BYT-ORB781999

Artikelname	General 2-Arachidonoylglycerol (2-AG) ELISA Kit
Artikelnummer	BYT-ORB781999
Hersteller Artikelnummer	orb781999
Alternativnummer	BYT-ORB781999-48, BYT-ORB781999-96
Hersteller	Biorbyt
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	All
Produktbeschreibung	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with 2-Arachidonoylglycerol(2-AG) protein. Standards or samples are added to the appropriate microtiter plate we...
Konzentration	300 ng/mL
Detektionsbereich	4.69-300 ng/mL
Sensitivitaet	1.41 ng/mL
Proben	Serum, plasma and other biological fluids

Anwendungsbeschreibung

Application Notes: standard: 300 ng/mL. Test principle: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with 2-AG protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to 2-AG. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of 2-AG in the samples is then determined by comparing the OD of the samples to the standard curve