

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Transferrin (TF, Beta-1 Metal-binding Globulin, Serotransferrin, Siderophilin) (Not for Export EU), Rabbit USB-T8199-04U

| | |
|--------------------------|---|
| Artikelname | Transferrin (TF, Beta-1 Metal-binding Globulin, Serotransferrin, Siderophilin) (Not for Export EU), Rabbit |
| Artikelnummer | USB-T8199-04U |
| Hersteller Artikelnummer | T8199-04U |
| Alternativnummer | USB-T8199-04U-1 |
| Hersteller | US Biological |
| Wirt | Rabbit |
| Kategorie | Antikörper |
| Immunogen | Purified transferrin isolated from pooled horse serum. |
| Produktbeschreibung | Applications: Suitable for use in IEP and RID. Other applications have not been tested. Recommended Dilutions: IEP: Use 2ul serum, plasma or equivalent against 120ul antiserum RID: Use a rosette arrangement with 10ul antiserum in 3 mm diameter centre... |
| Reinheit | Serum |
| Formulierung | Supplied as a lyophilized powder. No preservatives added. Reconstitute with 1ml sterile ddH ₂ O. Dilute further with PBS, pH 7.2. |