

Product Datasheet

Trident Membrane Protein Extraction Kit GTX-GTX16373

Artikelname	Trident Membrane Protein Extraction Kit
Artikelnummer	GTX-GTX16373
Hersteller Artikelnummer	GTX16373
Alternativnummer	GTX16373-5,GTX16373-20
Hersteller	Genetex
Kategorie	Kits/Assays
Applikation	ELISA, IP, WB
Anwendungsbeschreibung	1. Read the entire protocol carefully. 2. Before experiment, thaw buffer A or/and buffer B completely, invert the bottles a few times and keep them on ice. Place the filter cartridges in collection tubs, and pre-chill on ice. Pre-chill PBS on ice. <