

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### RM B-NGF SHD-200-11-1MG

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | RM B-NGF                                                                                                                                                                                                                                                      |
| Artikelnummer            | SHD-200-11-1MG                                                                                                                                                                                                                                                |
| Hersteller Artikelnummer | 200-11-1MG                                                                                                                                                                                                                                                    |
| Alternativnummer         | SHD-200-11-1MG                                                                                                                                                                                                                                                |
| Hersteller               | Fujifilm Irvine Scientific                                                                                                                                                                                                                                    |
| Kategorie                | Sonstiges                                                                                                                                                                                                                                                     |
| Applikation              | Bioassay, SDS-PAGE                                                                                                                                                                                                                                            |
| Spezies Reaktivität      | Mouse                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | Nerve growth factor beta (beta-NGF) is a neurotrophic factor that is important for the development and maintenance of sensory and sympathetic neurons. beta-NGF signals through the low affinity nerve growth factor receptor (LNGFR) and the tropomyosin ... |
| Molekulargewicht         | 13.6/27.2                                                                                                                                                                                                                                                     |
| UniProt                  | <a href="#">P01139</a>                                                                                                                                                                                                                                        |
| Puffer                   | 0.1% Trifluoroacetic Acid (TFA)                                                                                                                                                                                                                               |
| Quelle                   | E.coli                                                                                                                                                                                                                                                        |
| Reinheit                 | 95                                                                                                                                                                                                                                                            |
| Formulierung             | Lyophilized from a 0.2 µm filtered solution                                                                                                                                                                                                                   |

|                        |                                                                                                                                                                                                                                                                |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sequenz                | MSSTHPVFHM GEFSVCDSVS VWVGDKTTAT DIKGKEVTVL AEVNINNSVF<br>RQYFFETKCR ASNPVESGCR GIDSKHWNSY CTTHTFVKA<br>LTTDEKQAAW RFIRIDTACV CVLSRKATRR G                                                                                                                     |
| Anwendungsbeschreibung | Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0. |