

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

CoralHue humanized-codon dimeric Keima-Red (phdKeima-Red-MNLinker) MBL-AM-V0100M

Article Name	CoralHue humanized-codon dimeric Keima-Red (phdKeima-Red-MNLinker)
Biozol Catalog Number	MBL-AM-V0100M
Supplier Catalog Number	AM-V0100M
Alternative Catalog Number	MBL-AM-V0100M
Manufacturer	MBL
Category	Biochemikalien
Source	The CoralHue dKeima-Red gene was originally cloned from the stony coral (Montipora sp.).
Purity	A260/A280 > 1.5
Form	Dry form. Reconstitute with distilled water or TE before use.
Storage	-20°C
Formula	Dry form. Reconstitute with distilled water or TE before use.
Target	Keima-Red (phdKeima-Red-MNLinker)
Notes	For Research Use Only. Not for use in diagnostic procedures.