



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

3T3 cells, Unconjugated NMB-X1005

Article Name	3T3 cells, Unconjugated
Biozol Catalog Number	NMB-X1005
Supplier Catalog Number	X1005
Alternative Catalog Number	NMB-X1005
Manufacturer	NordicMubio
Category	Proteine/Peptide
Application	WB
Conjugation	Unconjugated
Product Description	Total cell lysates are useful as both positive and negative controls in immunoblotting. A431, a human epidermoid carcinoma cell line, and the EGF-stimulated A431 lysates are used as negative and positive controls, respectively, when studying the phos
Concentration	See vial for concentration
Buffer	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Formula	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Application Notes	Suggested use: A volume of 10-25ul is suggested for immunobloting.