

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

FIX&PERM Cell Fixation and Permeabilization Kit (CE/IVD) **NMB-GAS-002**

| | |
|----------------------------|--|
| Article Name | FIX&PERM Cell Fixation and Permeabilization Kit (CE/IVD) |
| Biozol Catalog Number | NMB-GAS-002 |
| Supplier Catalog Number | GAS-002 |
| Alternative Catalog Number | NMB-GAS-002 |
| Manufacturer | NordicMubio |
| Category | Sonstiges |
| Application | FC |
| Product Description | Flow cytometric analyses with monoclonal antibodies were so far mainly restricted to cell surface molecules. Intracellular structures such as cytoplasmic or nuclear enzymes, oncoproteins, cytokines, immunoglobulins etc. were largely excluded from suc... |
| Buffer | This FIX&PERM Cell Fixation and Permeabilization Kit contains 2 reagents: Fixation Medium (Reagent A) and Permeabilization Medium (Reagent B). It is intended for first fixing cells in suspension with Reagent A and then permeabilizing the cell membrane |
| Application Notes | Biological fluids (blood, bone marrow, and others) must be collected under sterile conditions. Anticoagulation with EDTA or heparin is recommended. The samples should be stored at room temperature until used. For optimal results, samples should be proces |