

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Goat IgG anti-Guinea Pig IgG2 (Fc)-FITC, MinX none NMB-GAGP/IGG2/FITC

Article Name	Goat IgG anti-Guinea Pig IgG2 (Fc)-FITC, MinX none
Biozol Catalog Number	NMB-GAGP/IGG2/FITC
Supplier Catalog Number	GAGp/lgG2/FITC
Alternative Catalog Number	NMB-GAGP/IGG2/FITC
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Species Reactivity	Guinea pig
Conjugation	FITC
Format	IgG
Target Specificity	IgG2 (Fc)
Cross-Adsorption (MinX)	no cross-adsorbtion
Product Description	The reactivity of the antiserum is directed to the subclass IgG2. It does not react with IgG1, IgG/Fab fragments, IgM and IgA or any non-Ig protein in guinea pig serum, as tested by immunoelectrophoresis and double radial immunodiffusion. To identify
Clonality	Polyclonal
Clone Designation	[Polyclonal]

Buffer	Fluorochrome-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water,
Source	Pools of highly purified homogenous IgG2 isolated from guinea pig serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Fluorochrome-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.
Antibody Type	Secondary Antibody
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence.