

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat IgG anti-Rabbit IgG+IgM+IgA (H+L)-unconj., MinX none NMB-GAR/IG/7S

Article Name	Goat IgG anti-Rabbit IgG+IgM+IgA (H+L)-unconj., MinX none
Biozol Catalog Number	NMB-GAR/IG/7S
Supplier Catalog Number	GAR/Ig/7S
Alternative Catalog Number	NMB-GAR/IG/7S
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Species Reactivity	Rabbit
Conjugation	Unconjugated
Format	IgG
Target Specificity	IgG+IgM+IgA (H+L)
Cross-Adsorption (MinX)	no cross-adsorbtion
Product Description	The reactivity of the antiserum is directed to the major isotypes of the rabbit immunoglobulin system (classes and both light chain types) including antibodies to common determinants, to class and to the surface determinants of the common Fab portion...
Clonality	Polyclonal
Clone Designation	[Polyclonal]

Buffer	Purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove i
Source	Highly purified IgG and pools of homogenous IgA and IgM isolated from rabbit serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formula	Purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.
Antibody Type	Secondary Antibody
Application Notes	Indirect immunofluorescence,ELISA,Dot blot,Immunoblotting.