

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat Peroxidase-anti-Peroxidase, Clone: [Polyclonal], Monoclonal NMB-G/PAP

Article Name	Goat Peroxidase-anti-Peroxidase, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-G/PAP
Supplier Catalog Number	G/PAP
Alternative Catalog Number	NMB-G/PAP
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Product Description	Peroxidase staining of unlabelled goat IgG antibodies
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Hyperimmune goat IgG against horseradish peroxidase reacted with peroxidase lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). 0.005% Merthiolate. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insolu
Source	Highly purified horseradish peroxidase.
Formula	Hyperimmune goat IgG against horseradish peroxidase reacted with peroxidase lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). 0.005% Merthiolate.