

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## NM LYSE: Flow Cytometry Lysing Solution (CE/IVD) NMB-GAS-003-1

Article Name	NM LYSE: Flow Cytometry Lysing Solution (CE/IVD)
Biozol Catalog Number	NMB-GAS-003-1
Supplier Catalog Number	GAS-003-1
Alternative Catalog Number	NMB-GAS-003-1
Manufacturer	NordicMubio
Category	Sonstiges
Application	FC
Product Description	For Wash- or No-Wash Lysing Procedures with Whole Blood or Marrow Samples NM-LYSE is a premixed, ready to use lysing solution fomulated for lysing erythrocytes following monoclonal antibody staining of whole blood. Treatment with this reagent simulta
Application Notes	No-Wash Staining and Lysing Procedure - For each sample add 50 $\mu$ l of EDTA anti-coagulated blood to a 3-5 ml tube - Add 20 $\mu$ l of the appropriate Nordic-MUbio monoclonal antibody conjugate - Incubate the tube for 15 minutes at 4C or at room temperature in

## BIOZOL