

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Rabbit anti Bacillus thermoproteolyticus rokko Thermolysin, Clone: [Polyclonal], Monoclonal NMB-NE153/PAB**

|                            |   |
|----------------------------|---|
| Article Name               | Rabbit anti Bacillus thermoproteolyticus rokko Thermolysin, Clone: [Polyclonal], Monoclonal   |
| Biozol Catalog Number      | NMB-NE153/PAB   |
| Supplier Catalog Number    | NE153/PAb   |
| Alternative Catalog Number | NMB-NE153/PAB   |
| Manufacturer               | NordicMubio   |
| Host                       | Rabbit  |
| Category                   | Antikörper  |
| Product Description        | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...   |
| Clonality                  | Monoclonal  |
| Clone Designation          | [Polyclonal]  |
| UniProt                    | <a href="#">P00800</a>  |
| Buffer                     | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitute |

|                   |  |
|-------------------|--|
| Source            | Thermolysin isolated and purified from <i>Bacillus thermoproteolyticus</i> rokko. Freund's complete adjuvant is used in the first step of the immunization procedure.  |
| Formula           | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Application Notes | Immunoprecipitation.   |