

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rabbit IgG anti-Sheep IgG (H+L)-HRPO, MinX none NMB-RASH/IGG(H+L)/PO

|                            |   |
|----------------------------|---|
| Article Name               | Rabbit IgG anti-Sheep IgG (H+L)-HRPO, MinX none   |
| Biozol Catalog Number      | NMB-RASH/IGG(H+L)/PO  |
| Supplier Catalog Number    | RASH/IgG(H+L)/PO  |
| Alternative Catalog Number | NMB-RASH/IGG(H+L)/PO  |
| Manufacturer               | NordicMubio   |
| Host                       | Rabbit  |
| Category                   | Antikörper  |
| Species Reactivity         | Sheep   |
| Conjugation                | HRPO  |
| Format                     | IgG   |
| Target Specificity         | IgG (H+L)   |
| Cross-Adsorption (MinX)    | no cross-adsorbtion   |
| Product Description        | The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in sheep serum, as... |
| Clonality                  | Polyclonal  |
| Clone Designation          | [Polyclonal]  |

|                   |   |
|-------------------|---|
| Buffer            | Peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, |
| Source            | Purified normal IgG isolated from pooled sheep serum. Freund's complete adjuvant is used in the first step of the immunization procedure.   |
| Form              | HRP   |
| Formula           | Peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.   |
| Antibody Type     | Secondary Antibody  |
| Application Notes | ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.  |