



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## 8-iso-PGF2alpha ELISA kit ENZ-ADI-900-010

A.C. I. Ni.	0.1. 0050 1.1. 5154.13
Article Name	8-iso-PGF2alpha ELISA kit
Biozol Catalog Number	ENZ-ADI-900-010
Supplier Catalog Number	ADI-900-010
Alternative Catalog Number	ENZ-ADI-900-010-96
Manufacturer	Enzo Life Sciences
Category	Kits/Assays
Application	ELISA
Species Reactivity	All
Product Description	The measurement of F2-isoprostanes has emerged as one of the most reliable approaches to assess oxidative stress status in vivo, providing an important tool to explore the role of oxidative stress in the pathogenesis of human disease. This is a color
Sensitivity	16.3 pg/ml (range 6.1 - 100,000 pg/ml)
Target	8-iso-PGF2 alpha
Application Notes	For the quantitative determination of 8-iso-PGF2alpha in culture supernatants, tissue, and urine from any species. Cited sample types include plasma.