

Diagnostica Vertrieb GmbH, Leipziger Straße 4 85386 Eching, Germany **Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99 **E-Mail:** info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## DNA Marker Low ABN-D0001

| Article Name               | DNA Marker Low   |
|----------------------------|--|
| Biozol Catalog Number      | ABN-D0001  |
| Supplier Catalog Number    | D0001  |
| Alternative Catalog Number | ABN-D0001-1  |
| Manufacturer               | Abnova   |
| Category                   | Molekularbiologie  |
| Application                | EP, PCR  |
| Product Description        | DNA Marker Low is useful for the determination of DNA fragments 50 bp to 1,000 bp in size. The DNA size marker consists of 10 double-stranded, blunt-end DNA fragments. The DNA bands close to each other, 550 bp and 500 bp, 120 bp and 100 bp, are for r |
| Application Dilute         | DNA sample: 5 uL6 X BPB loading dye: 1 uL (addition of 1/5 volume of DNA sample)Mix and load to a well   |