



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

DNA Marker High ABN-D0002

Article Name	DNA Marker High
Biozol Catalog Number	ABN-D0002
Supplier Catalog Number	D0002
Alternative Catalog Number	ABN-D0002-1
Manufacturer	Abnova
Category	Molekularbiologie
Application	EP, PCR
Product Description	DNA Marker High is useful for the determination of DNA fragments 300 bp to 10,000 bp in size. The DNA size marker consists of 13 double-stranded, blunt-end DNA fragments. The DNA bands close to each other, 525 bp and 500 bp, are for reference indicat
Application Dilute	DNA sample: 5 uL6 X BPB loading dye: 1 uL (addition of 1/5 volume of DNA sample)Mix and load to a well