



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Fastback Epoxy Activated 4 High Flow Resin (50ml) PAR-FASTBACK-EPOXY-4HF-50

Article Name	Fastback Epoxy Activated 4 High Flow Resin (50ml)
Biozol Catalog Number	PAR-FASTBACK-EPOXY-4HF-50
Supplier Catalog Number	Fastback-EPOXY-4HF-50
Alternative Catalog Number	PAR-FASTBACK-EPOXY-4HF-50
Manufacturer	ProteinArk
Category	Sonstiges
Product Description	Fastback Epoxy Activated 4 High Flow Resin (50ml)