

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Trident Membrane Protein Extraction Kit GTX16373

Article Name	Trident Membrane Protein Extraction Kit
Biozol Catalog Number	GTX16373
Supplier Catalog Number	GTX16373
Alternative Catalog Number	GTX16373-20,GTX16373-5
Manufacturer	GeneTex
Category	Sonstiges
Application	ELISA, IP, WB
Application Notes	1. Read the entire protocol carefully. 2. Before experiment, thaw buffer A or/and buffer B completely, invert the bottles a few times and keep them on ice. Place the filter cartridges in collection tubs, and pre-chill on ice. Pre-chill PBS on ice.