



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Fixation/Permeabilization, (buffer) ENZ-ADI-950-010-100

Article Name	Fixation/Permeabilization, (buffer)
Biozol Catalog Number	ENZ-ADI-950-010-100
Supplier Catalog Number	ADI-950-010-100
Alternative Catalog Number	ENZ-ADI-950-010-100-100
Manufacturer	Enzo Life Sciences
Category	Sonstiges
Application	FC
Product Description	Fixation/Permeabilization, (buffer)
Purity	Purified
Form	Liquid. PBS, pH 7.2, containing 4% formaldehyde and 0.1% saponin.
Application Notes	Intended for use in simultaneously fixing and permeabilizing cells in one step to facilitate immunostaining and subsequent analysis of intracellular proteins using flow cytometry. This buffer can be used as an alternative to a two step process with fixat