



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Permeabilization, (buffer) ENZ-ADI-950-012-100

Article Name	Permeabilization, (buffer)
Biozol Catalog Number	ENZ-ADI-950-012-100
Supplier Catalog Number	ADI-950-012-100
Alternative Catalog Number	ENZ-ADI-950-012-100-100
Manufacturer	Enzo Life Sciences
Category	Sonstiges
Application	FC
Product Description	Permeabilization, (buffer)
Purity	Purified
Form	Liquid. PBS, pH 7.2, containing 0.1% saponin and 0.09% sodium azide.
Application Notes	Intended for use in permeabilizing cells to facilitate immunostaining and subsequent analysis of intracellular proteins using flow cytometry. Use of this buffer should be preceded by fixation of the cells using a suitable fixation buffer (Prod. No. ADI-9