

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

pro-SUMO-1 (human), (recombinant) (His-tag) ENZ-BML-UW9190-0500

Article Name	pro-SUMO-1 (human), (recombinant) (His-tag)
Biozol Catalog Number	ENZ-BML-UW9190-0500
Supplier Catalog Number	BML-UW9190-0500
Alternative Catalog Number	ENZ-BML-UW9190-0500-500
Manufacturer	Enzo Life Sciences
Category	Proteine/Peptide
Product Description	pro-SUMO-1 (human), (recombinant) (His-tag)...
Molecular Weight	16.4 kDa
UniProt	P63165
Source	Produced in E. coli. Human pro-SUMO-1 is fused at the N-terminus to a His-tag.
Purity	95% (SDS-PAGE)
Form	Liquid. In 20mM HEPES, pH 8.0, 50mM NaCl, 1mM DTT.
Target	SUMO
Application Notes	Uses: 1. Study of SUMO-1 regulation and processing.2. Substrate for deSUMOylating enzymes (C-terminal hydrolase activity).