

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **SUMO-3 (human), (recombinant) (His-tag) ENZ-BML-UW9215-0500**

Article Name	SUMO-3 (human), (recombinant) (His-tag)
Biozol Catalog Number	ENZ-BML-UW9215-0500
Supplier Catalog Number	BML-UW9215-0500
Alternative Catalog Number	ENZ-BML-UW9215-0500-500
Manufacturer	Enzo Life Sciences
Category	Proteine/Peptide
Product Description	SUMO-3 (human), (recombinant) (His-tag)...
Molecular Weight	15.5 kDa
UniProt	<a href="#">P55854</a>
Source	Produced in E. coli. Human SUMO-3 is fused to an N-terminal His-tag.
Purity	95% (SDS-PAGE)
Form	Liquid. In 20mM HEPES, pH 8.0, containing 50mM sodium chloride and 1mM dithiothreitol.
Target	SUMO
Application Notes	Uses: 1. SUMO-3 modification of specific proteins in vitro. 2. Demonstrate novel proteins are potential targets for SUMOylation under in vitro conditions. 3. Generate substrates for deSUMOylating enzymes.