



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## AMPIPROBE UP/Mg/Mh Extraction Controls ENZ-ENZ-GEN210-0500

Article Name	AMPIPROBE UP/Mg/Mh Extraction Controls
	/ with the self-order of the s
Biozol Catalog Number	ENZ-ENZ-GEN210-0500
Supplier Catalog Number	ENZ-GEN210-0500
Alternative Catalog Number	ENZ-ENZ-GEN210-0500-500
Manufacturer	Enzo Life Sciences
Category	Sonstiges
Species Reactivity	Bacteria
Product Description	AMPIPROBE UP/Mg/Mh Extraction Controls
Target	Ureaplasma spp.  Mycoplasma genitalium  Mycoplasma hominis
Application Notes	The AMPIPROBE UP/Mg/Mh Extraction Controls is compatible with manual and automated DNA extraction methods as well as any properly calibrated qPCR thermal cycler capable to detect fluorescence decay. It has been validated for use on the QIAGEN QIAsymphon