

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Triple FLAG Peptide ROC-000-000-H47

|                            |   |
|----------------------------|---|
| Article Name               | Triple FLAG Peptide   |
| Biozol Catalog Number      | ROC-000-000-H47   |
| Supplier Catalog Number    | 000-000-H47   |
| Alternative Catalog Number | ROC-000-000-H47   |
| Manufacturer               | Rockland Immunochemicals  |
| Category                   | Proteine/Peptide  |
| Application                | SDS-PAGE  |
| Conjugation                | Unconjugated  |
| Concentration              | 1.0 mg/ml by dry weight   |
| Buffer                     | None  |
| Form                       | Lyophilized   |
| Application Dilute         | IP: YES - see PROTOCOL, WB: YES as blocking peptide in PCA  |
| Application Notes          | The Triple FLAG Peptide has been tested by SDS-PAGE. A working concentration of 500 µg/mL for triple FLAG peptide is typically used for competitive elution of FLAG containing proteins from anti-DYKDDDDK affinity gels (p/n 200-350-383). For use as a co |