



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Collagen I for plate coating - Bovine Skin ROC-001-G70-MF6

Article Name	Collagen I for plate coating - Bovine Skin
Biozol Catalog Number	ROC-001-G70-MF6
Supplier Catalog Number	001-G70-MF6
Alternative Catalog Number	ROC-001-G70-MF6
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Application	SDS-PAGE
Conjugation	Unconjugated
Concentration	1.0 mg/mL by dry weight
Buffer	0.013M HCI
Form	Liquid (sterile filtered)
Application Notes	Bovine Collagen Type I solution has been tested in SDS-page and is suitable for plate coating. Enough to coat 10 (ten) T75 plates.