

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit IgG anti-Dog IgG (F(ab)2)-FITC, MinX none ROC-204-4204

| | |
|----------------------------|---|
| Article Name | Rabbit IgG anti-Dog IgG (F(ab)2)-FITC, MinX none |
| Biozol Catalog Number | ROC-204-4204 |
| Supplier Catalog Number | 204-4204 |
| Alternative Catalog Number | ROC-204-4204 |
| Manufacturer | Rockland Immunochemicals |
| Host | Rabbit |
| Category | Antikörper |
| Application | FLISA,FACS,IF |
| Species Reactivity | Canine |
| Immunogen | Dog IgG F(ab)2 fragment |
| Conjugation | FITC |
| Format | IgG |
| Target Specificity | IgG (F(ab')2) |
| Cross-Adsorption (MinX) | no cross-adsorbtion |
| Clonality | Polyclonal |
| Concentration | 10.0 mg/mL by UV absorbance at 280 nm |
| Buffer | 0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Form | Lyophilized |

| | |
|--------------------|---|
| Target | Dog |
| Antibody Type | Secondary Antibody |
| Application Dilute | FLISA: 1:10,000 - 1:50,000, Flow Cytometry: 1:500 - 1:2,500, IF Microscopy: 1:1,000 - 1:5,000 |
| Application Notes | Secondary antibody reagents are ideal for ELISA, western blotting, Immunohistochemistry, Fluorescence Microscopy, Flow Cytometry as well as other antibody detection methods. |